

MICRONIC

Traceable Sampling Solutions

Micronic Expell Plus™ Tips

Manufactured from special polymers and with an extreme smooth inside, the Micronic Expell Plus™ Tips have been developed for a wide range of applications in the field of liquid handling. The use of high-grade pure polypropylene virtually eliminates samples retention and reduces the binding of DNA, RNA and proteins to a very minimum.

The comprehensive line of tips are fully compatible with the Micronic (Electronic) Pipettors, they provides highest precision and accuracy without any high valuable sample loss.



Standard Tip



ExpellPlus™ low retention

Plus™ low retention properties

Tests show that the Expell Plus™ tips deliver volumes within 0,1% of the set volume, versus 0,5-0,7% for standard tips. This significant difference in sample retention also is clearly visible when pipetting a dyed liquid (see left)

The Expell™ tips without the Plus™ low retention properties are of the same high manufacturing standards, maintaining strict tolerances with ideal wetting properties and high transparency and purity levels.

Micronic Expell Plus™ Filter tips

The Micronic Expell Plus™ Filter tips are specially designed to prevent airborne contaminants from coming in contact with the barrel of the pipettor. The Micronic Expell Plus™ Filter tips are tested and certified to conform to strict bio-load specifications:



- No detectable nucleic acid contamination
- No detectable endonucleases
- No detectable endotoxin
- Sterile




Easily accessible, space and waste saving

The racked Expell™ tips are supplied in handy dispenser boxes occupying minimal space in the bench, with easy access to the next rack to be used. Due to innovative packaging solutions, the racked tips take up less storage space than conventional tips and are waste saving at the same time

For more information about our products contact NovaVeth B.V.

NovaVeth B.V.
Huizermaatweg 344
1276 LJ Huizen
The Netherlands

T: +31 (0)35 525 88 82
F: +31 (0)35 525 81 55
E: sales@novaveth.com
W: www.novaveth.com









MICRONIC

Traceable Sampling Solutions






Micronic *Expell Plus™* Tips

Order information

Expell™ and Expell Plus™ Filter tips

	Cat#:	Description:	Qty:	Packaging:
	E5030030	Expell Plus™ 10µl FT	10 racks of 96	Racked/Sterile
	E5030060	Expell Plus™ 10µl FT, long	10 racks of 96	Racked/Sterile
	E5030062	Expell Plus™ 20µl FT	10 racks of 96	Racked/Sterile
	E5030066	Expell Plus™ 100µl FT	10 racks of 96	Racked/Sterile
	E5030090	Expell Plus™ 200µl FT	10 racks of 96	Racked/Sterile
	E5130150	Expell™ 1000µl FT	8 racks of 60	Racked/Sterile

Expell™ and Expell Plus™ Pipette tips

	Cat#:	Description:	Qty:	Packaging:
	E5030010	Expell Plus™ 10µl	1000	Bulk
	E5030015	Expell Plus™ 10µl	10 racks of 96	Racked
	E5030020	Expell Plus™ 10µl	10 racks of 96	Racked/Sterile
	E5130040	Expell Plus™ 10µl, long	1000	Bulk
	E5130045	Expell Plus™ 10µl, long	10 racks of 96	Racked
	E5130050	Expell Plus™ 10µl, long	10 racks of 96	Racked/Sterile
	E5130070	Expell Plus™ 200µl	1000	Bulk
	E5130075	Expell Plus™ 200µl	10 racks of 96	Racked
	E5130080	Expell Plus™ 200µl	10 racks of 96	Racked/Sterile
	E5130100	Expell™ 300µl	1000	Bulk
	E5130105	Expell™ 300µl	10 racks of 96	Racked
	E5130110	Expell™ 300µl	10 racks of 96	Racked/Sterile
	E5130130	Expell™ 1000µl	1000	Bulk
	E5130135	Expell™ 1000µl	16 racks of 60	Racked
	E5130140	Expell™ 1000µl	16 racks of 60	Racked/Sterile

For more information about our products contact NovaVeth B.V.

NovaVeth B.V.
 Huizermaatweg 344
 1276 LJ Huizen
 The Netherlands

T: +31 (0)35 525 88 82
 F: +31 (0)35 525 81 55
 E: sales@novaveth.com
 W: www.novaveth.com